


CRYO-CYL

LIQUID CYLINDER

Like the Dura-Cyl models, the Cryo-Cyl 80HP model is designed and built to meet the rugged demands of the liquid cylinder market. However, in contrast, this model is designed specifically for liquid and low to medium gas flow applications. By specifically targeting these applications, we are able to offer this model at an economical value over our premium Dura-Cyl series.



PRODUCT HIGHLIGHTS

- Ideal for liquid nitrogen, oxygen, argon, CO₂ or nitrous oxide
- Stainless steel construction
- Thick, dent-resistant outer shell
- Patented durable, inner-vessel support system
- Heavy-duty footing and large diameter handling ring with two supports
- Roto-Tel Liquid Level Gauge System
- Five-year vacuum warranty
- Vent and Gas use valve are combined for simple and safe operation
- Built-in economizer circuit prevents gas loss during gas intermittent use




Innovation. Experience. Performance.®

CRYO-CYL

LIQUID CYLINDER

Specifications		
MODEL		80 L
	Pressure	HP
	Part Number	10648610
CAPACITY⁽¹⁾⁽²⁾		
Liquid (Gross)	(liters)	85
Liquid (Net)	(liters)	80
Gas (N ₂)	ft ³ / Nm ³	1,680 / 44
Gas (O ₂)	ft ³ / Nm ³	2,108 / 55
Gas (Ar)	ft ³ / Nm ³	2,049 / 54
Gas (CO ₂)	ft ³ / Nm ³	1,640 / 43
Gas (N ₂ O)	ft ³ / Nm ³	1,555 / 41
PERFORMANCE		
NER (N ₂)	% per day	3.0
NER (O ₂ - Ar)	% per day	2.0
NER (CO ₂ - N ₂ O)	% per day	0.8
Gas Flow (N ₂ , O ₂ , Ar) ⁽³⁾	SCFH/Nm ³ /hr	100 / 2.6
Gas Flow (CO ₂ , N ₂ O)	SCFH/Nm ³ /hr	35 / 0.9
DIMENSIONS & PRESSURE RATINGS		
Relief Valve Setting	psig / barg	350 / 24
Operating Pressure ⁽⁴⁾	psig / barg	125 / 8.6
DOT/TC Rating		4L292
Diameter	in / cm	20 / 50.8
Height ⁽⁵⁾	in / cm	39.5 / 100.3
Tare Weight	lb / kg	165 / 74.8
Full Weight (N ₂)	lb / kg	287 / 130
(O ₂)	lb / kg	340 / 155
(Ar)	lb / kg	377 / 171
(CO ₂)	lb / kg	353 / 161



Nomenclature

1. Gas Use / Vent Valve
2. Fill / Liquid Valve
3. Pressure Building Valve
4. Pressure Building Regulator
5. Pressure Gauge
6. Pressure Relief Valve
7. Rupture Disk
8. Liquid Level Gauge

Model: HP - High Pressure (Liquid or Gas Use)

- (1) Net gas capacities at DOT 4L limits.
- (2) The Cryo-Cyl model is available with permanently installed CGA fittings for medical applications. Contact Customer Service for details.
- (3) Gas flows of twice the continuous flow rate can be achieved for one hr. over an eight hr. period.
- (4) Pressure building regulator range (50-175 psi).
- (5) Height dimensions are measured from the floor to the top of the sight gauge protector.

Chart Inc.
U.S.: 1-800-400-4683
Worldwide: 1-952-758-4484



© 2013 Chart Inc.
P/N 11035057
www.chartindustries.com